

### Product Specification Sheet

Description: CypExpress™ 2B6  
Product Number: CE2B6.250  
Aliquot: 250 mg  
Lot Number: ce2B6.160915  
Storage: -80°C

- DESCRIPTION:** CypExpress™ 2B6 is a permeabilized and stabilized dried yeast powder preparation containing recombinant human CYP 2B6 and recombinant human P450 NADPH oxidoreductase.
- ACTIVITY:** A reaction mixture containing 20 µM 7-ethoxy-4-(trifluoromethyl)coumarin and 100 mg/mL of CypExpress™ 2B6 produced 2.6 uM 7-Hydroxy-4-(trifluoromethyl)coumarin in 3 hours at 30°C under the conditions described.
- THAWING:** It is important to avoid water adsorption during thawing. Remove the CypExpress™ vial from storage and warm to room temperature prior to opening. Return unused powder to -80°C.
- STORAGE:** Store at -80°C. Although CypExpress™ preparations can be warmed to room temperature multiple times, it is best to aliquot as needed to avoid multiple freeze/thaw cycles.
- SUGGESTED PROTOCOLS:** The basic protocols attached are intended as a starting point that can be optimized by an experienced user for specific applications and samples.