

**Product Specification Sheet**

Description: CypExpress™ 2D6  
Product Number: CE2D6.10  
Aliquot: 10 Grams  
Lot Number: ce2d6.2  
Storage: -80°C

- DESCRIPTION:** CypExpress™ 2D6 is a permeabilized and stabilized dried yeast powder preparation containing recombinant human CYP 2D6 and recombinant human P450 NADPH oxidoreductase.
- ACTIVITY:** A reaction mixture containing 1 mM dextromethorphan and 100 mg/mL of CypExpress™ 2D6 will produce 268.8 uM dextrophan in 4 hours at 30°C under the conditions described.
- THAWING:** It is important to avoid water adsorption during thawing. Remove the CypExpress™ vial from storage and warm to room temperature prior to opening. Return unused powder to -80°C.
- STORAGE:** Store at -80°C. Although CypExpress™ preparations can be warmed to room temperature multiple times, it is best to aliquot as needed to avoid multiple freeze/thaw cycles.
- SUGGESTED PROTOCOLS:** The basic protocols attached are intended as a starting point that can be optimized by an experienced user for specific applications and samples.