

**Product Specification Sheet**

Description: Human Factor Va  
Product Number: CF23  
Aliquot: 0.05 mg  
Lot Number: cf23.3340  
Storage: -20°C

<b>DESCRIPTION:</b>	Human Factor Va derived from factor V treated with thrombin and immunoaffinity purified. This prep is free of activation peptides and thrombin.
<b>MOLECULAR WEIGHT:</b>	168 kDa
<b>CONCENTRATION:</b>	3.6 mg/ml
<b>STORAGE BUFFER:</b>	50% Glycerol
<b>STORAGE:</b>	Store at -20°C or colder.
<b>SPECIFIC ACTIVITY:</b>	2120 U/mg
<b>SOURCE:</b>	Human plasma
<b>PURIFICATION:</b>	Factor Va is prepared by activating purified factor V with thrombin and then purified by immunoaffinity chromatography.
<b>PURITY:</b>	>95% by SDS PAGE Analysis