

Product Specification Sheet

Description: Anti-human SCF
monoclonal antibody
Product Number: CK51
Aliquot: 1.0 mg
Lot Number: Typical Lot
Storage: -80°C

DESCRIPTION:	Monoclonal antibody against Human Stem Cell Factor (SCF)
MOLECULAR WEIGHT:	18.5 kDa
STORAGE BUFFER:	0.5 x PBS, pH 7.2
STORAGE:	Once reconstituted, aliquot and store at -80°C.
IMMUNOGEN:	Recombinant Human SCF
PURIFICATION:	IgG1 antibody was purified by Protein G affinity chromatography.
HOST:	Mouse
CLASS:	IgG1
APPLICATIONS:	ELISA: 4.0 to 8.0 ug/mL Western: 0.25 to 0.5 ug/mL
SPECIES SPECIFICITY:	Human and Rat