

Product Specification Sheet

Description: Anti-human BMP-7
monoclonal antibody
Product Number: CK67
Aliquot: 1.0 mg
Lot Number: Typical Lot
Storage: -80°C

DESCRIPTION: Monoclonal antibody against Human Bone Morphogenetic Protein 7 (BMP-7)

MOLECULAR WEIGHT: 29 kDa

STORAGE BUFFER: 0.5 x PBS, pH 7.2

STORAGE: Once reconstituted, aliquot and store at -80°C.

IMMUNOGEN: Recombinant Human BMP-7

PURIFICATION: IgG1 antibody was purified by Protein G affinity chromatography.

HOST: Mouse

CLASS: IgG1

APPLICATIONS: ELISA: 1.0 ug/mL
Western: 2.0 ug/mL

SPECIES SPECIFICITY: Human