



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Bovine Fibrinogen
Product Number: CT28
Aliquot: 25 mg
Lot Number: Typical Lot
Storage: -70°C

DESCRIPTION:	Prepared from fresh bovine plasma using several chromatographic steps.
MOLECULAR WEIGHT:	330,000
CONCENTRATION:	40 mg/mL
SIZE:	25.0 mg
VOLUME:	0.625 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Bovine plasma
BUFFER:	0.02M Sodium Citrate-HCl; pH 7.4

For research purposes only.