



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Dog Fibrinogen
Product Number: CT29
Aliquot: 1.0 mg
Lot Number: Typical Lot
Storage: -70°C

DESCRIPTION:	Prepared from fresh dog plasma using several chromatographic steps.
MOLECULAR WEIGHT:	330,000
CONCENTRATION:	40 mg/mL
SIZE:	1.0 mg
VOLUME:	0.025 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Dog plasma
BUFFER:	0.02M Sodium Citrate-HCl; pH 7.4

For research purposes only.