



OXFORD BIOMEDICAL RESEARCH

## Product Specification Sheet

Rat Fibrinogen  
Product Number: CT31  
Aliquot: 1.0 mg  
Lot Number: Typical Lot  
Storage: -70°C

<b>DESCRIPTION:</b>	Prepared from fresh rat plasma using several chromatographic steps.
<b>MOLECULAR WEIGHT:</b>	340,000
<b>CONCENTRATION:</b>	10 mg/mL
<b>SIZE:</b>	1.0 mg
<b>VOLUME:</b>	0.1 mL
<b>STORAGE:</b>	-70°C
<b>FORM:</b>	Frozen liquid
<b>SOURCE:</b>	Rat plasma
<b>BUFFER:</b>	0.01M HEPES; 0.3M NaCl; 10 u/ml Trasylol; pH 7.4

For research purposes only.