

Product Specification Sheet

Anti-Human Estrogen Receptor β (ER β)
Product Number: ER01
Aliquot: 0.1 mg
Lot Number: er01.typical lot
Storage: -20° C

IMMUNOGLOBULIN:	Polyclonal Chicken IgY
IMMUNOGEN:	Recombinant Human ER β
PREDICTED SIZE:	53 kDa
CONCENTRATION:	1.0 mg/mL
SPECIFICITY:	Human ER β , not tested on other species.
STORAGE:	Store at -20° C. Avoid repeating freeze/thaw cycles, as it may result in loss of activity.
PURITY:	IgY
WORKING DILUTION:	Suggested starting concentration for Western blots is 1.0 μ g/mL. Not tested for IHC or ELISA.




Figure 1: Western blot of ER β using ECL detection. Final concentration of antibody was 1.0 μ g/mL.