

Product Specification Sheet

Anti-Human Estrogen Receptor β (ER β)
Product Number: ER02
Aliquot: 100 μ L
Lot Number: er02.typical lot
Storage: -20 $^{\circ}$ C

| | |
|--------------------------|--|
| IMMUNOGLOBULIN: | Polyclonal Rabbit IgG |
| IMMUNOGEN: | Recombinant Human ER β |
| PREDICTED SIZE: | 53 kDa |
| CONCENTRATION: | Undetermined |
| SPECIFICITY: | Human ER β , not tested on other species. |
| STORAGE: | Store at -20 $^{\circ}$ C. Avoid repeating freeze/thaw cycles, as it may result in loss of activity. |
| PURITY: | Serum |
| WORKING DILUTION: | Suggested working dilution for Western blots is 1:5000. Not tested for IHC or ELISA. |




Figure 1: Western blot of ER β using ECL detection. A dilution factor of 1:5000 was used.