

**Product Specification Sheet**

Mouse Fibrinogen, FITC Labeled  
Product Number: FB01  
Aliquot: 1.0 mg  
Lot Number: fb01.810  
Storage: -70°C

**DESCRIPTION:** Mouse fibrinogen that has been amino terminal labeled with FITC. Fibrinogen is prepared from fresh mouse plasma that has been plasminogen depleted.

**MOLECULAR WEIGHT:** 340,000 daltons

**CONCENTRATION:** 10.34 mg/mL

**STORAGE BUFFER:** 0.02M Sodium Citrate-HCl; pH 7.4

**STORAGE:** Store at -70°C and protect from light

**LABEL:** Labeled with FITC (fluorescein isothiocyanate) with an 11:1 mole ratio (FITC:Fibrinogen)

**PURIFICATION:** >95% pure by SDS-PAGE analysis

**THAWING:** Fibrinogen should be thawed at 37°C without agitation. Once it is completely thawed, GENTLY mix before use. Fibrinogen should be kept at 25-37°C while using. Flash freeze aliquots.