

Product Specification Sheet

Porcine Fibrinogen
Product Number: FB03
Aliquot: 25 mg
Lot Number: Typical Lot
Storage: -70°C

DESCRIPTION:	Porcine fibrinogen prepared from fresh plasma using chromatography. Plasminogen has been depleted by affinity chromatography.
MOLECULAR WEIGHT:	330 kDa
CONCENTRATION:	10 mg/mL
STORAGE BUFFER:	0.02M Sodium Citrate, pH 7.4
STORAGE:	-70°C
PURITY:	>95% by SDS-PAGE analysis
USAGE:	Thaw fibrinogen in a 37°C water bath. Aliquot as desired and freeze until needed. Keep at 25-37°C while using.