



OXFORD BIOMEDICAL RESEARCH

## Product Specification Sheet

Human Multimeric Vitronectin

Product Number: GP 02

Aliquot: 100  $\mu$ g

Lot Number: gp02.605A

Storage: -80°C

<b>DESCRIPTION:</b>	This human multimeric vitronectin is urea prepared from human plasma using non-denaturing chromatography.
<b>MOLECULAR WEIGHT:</b>	75 kDa
<b>CONCENTRATION:</b>	3.4 mg/mL determined spectrophotometrically.
<b>STORAGE BUFFER:</b>	0.05 M Sodium Phosphate; 0.1 M NaCl; pH 7.4
<b>STORAGE:</b>	-80°C. <b>AVOID REPEATED FREEZE THAWS.</b>