



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Human Vitronectin Antiserum

Product Number: GP 10

Aliquot: 1.0 mL

Lot Number: gp10.105A

Storage: -70°C

DESCRIPTION:	This rabbit polyclonal antibody is specific for human vitronectin and is provided as neat antiserum.
MOLECULAR WEIGHT:	75 kDa Monmeric 100 kDa Multimeric or Urea Prepared form
CONCENTRATION:	Not determined
STORAGE BUFFER:	Not applicable
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Vitronectin purified from human plasma using non-denaturing chromatography.
PURIFICATION:	Not purified
APPLICATIONS:	Suitable for use in Western blotting and ELISA applications.
SPECIES SPECIFICITY:	Human