



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Mouse Vitronectin Antiserum

Product Number: GP 14

Aliquot: 1.0 mL

Lot Number: gp14.1101A

Storage: -70°C

DESCRIPTION:	This rabbit polyclonal antibody is specific for mouse vitronectin and is provided as neat antiserum.
MOLECULAR WEIGHT:	75 kDa Monmeric 100 kDa Multimeric or Urea Prepared form
CONCENTRATION:	Not determined
STORAGE BUFFER:	Not applicable
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Vitronectin purified from mouse plasma using non-denaturing chromatography.
PURIFICATION:	Not purified
APPLICATIONS:	Suitable for use in Western blotting and ELISA applications. Also tested in immunohistochemistry 1:50 recommended dilution.
SPECIES SPECIFICITY:	Mouse