



## Product Specification Sheet

Rabbit Anti-Mouse Vitronectin  
(IgG Fraction)  
Product Number: GP16  
Aliquot: 1.0 mg  
Lot Number: gp16.105  
Storage: -80°C

<b>DESCRIPTION:</b>	Rabbit polyclonal antibody against vitronectin.
<b>MOLECULAR WEIGHT:</b>	75 kDa
<b>CONCENTRATION:</b>	3.4 mg/mL
<b>STORAGE BUFFER:</b>	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
<b>STORAGE:</b>	-80°C
<b>IMMUNOGEN:</b>	Vitronectin
<b>PURIFICATION:</b>	The IgG fraction was purified by immobilized Protein A.
<b>APPLICATIONS:</b>	ELISA, IHC, Western Blot
<b>SPECIES SPECIFICITY:</b>	Specific for mouse.