



OXFORD BIOMEDICAL RESEARCH

## Product Specification Sheet

Biotin Labeled Monomeric Human  
Vitronectin  
Product Number: GP31  
Aliquot: 0.1 mg  
Lot Number: Typical Lot  
Storage: -70°C

<b>DESCRIPTION:</b>	Prepared from fresh human plasma using non-denaturing chromatography, then biotin labeled at primary amines.
<b>MOLECULAR WEIGHT:</b>	75,000
<b>CONCENTRATION:</b>	1.0 mg/mL
<b>SIZE:</b>	0.1 mg
<b>VOLUME:</b>	0.1 mL
<b>STORAGE:</b>	-70°C
<b>FORM:</b>	Frozen liquid
<b>SOURCE:</b>	Human plasma
<b>BUFFER:</b>	0.05M Sodium Phosphate; 0.1M NaCl; pH 7.4

For research purposes only.