



Product Specification Sheet

Human Fibronectin
Product Number: GP 40
Aliquot: 1.0 mg
Lot Number gp40.406
Storage: -80°C

DESCRIPTION:	Human Fibronectin
CONCENTRATION:	22.0 mg/mL
MOLECULAR WEIGHT:	220,000 (reduced)
STORAGE BUFFER:	0.1M Tris-Hcl; 0.15M NaCl; pH 7.4
STORAGE:	Store at -80°C. Thaw at 37°C without agitation
SOURCE:	Prepared from fresh human plasma.
USES:	Ideal for tissue culture studies and protein-protein interactions.