



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Mouse Fibronectin
Product Number: GP41
Aliquot: 1.0 mg
Lot Number: Typical Lot
Storage: -70°C

DESCRIPTION:	Prepared from fresh mouse plasma. Ideal reagent for tissue culture studies and protein-protein interactions.
MOLECULAR WEIGHT:	220,000
CONCENTRATION:	10 mg/mL
SIZE:	1.0 mg
VOLUME:	0.1 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Mouse plasma
BUFFER:	0.1 M Tris-HCl; 0.15 M NaCl; pH 7.4

For research purposes only.