



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Rat GST Ya Antiserum

Product Number: GS 09

Aliquot: 100 μ L

Lot Number 071996DM

Storage: -70°C

DESCRIPTION:	Goat antisera specific for Rat Glutathione S Transferase Ya. This antibody does not cross react with the mu or pi class of GST.
CONCENTRATION:	Not Determined
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
STORAGE BUFFER:	None
IMMUNOGEN:	Affinity Purified Rat Liver GST Ya
PURIFICATION:	None
APPLICATIONS:	Western Blotting: 1:1000 This aniserum yields a positive result with 0.1 μ g of Purified Rat GST Ya when using a 1:1000 dilution of GS 09 and a 1:5000 dilution of alkaline phosphatase-conjugated anti-goat secondary antibody.
SPECIES SPECIFICITY:	Rat and Human.
REFERENCE:	Holmes, W., <i>et. al.</i> , <i>Br. J. Cancer</i> , 62 :209-212 (1990)