



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA
USA: 800-692-4633 • Fax: 248-852-4466
www.oxfordbiomed.com

Anti-Rat GST Yb

Product # GS 23

Lot # 112596DM

- Contents:** Polyclonal goat antiserum
- Immunogen:** Affinity purified rat GST Yb1.
- Specificity:** WESTERN BLOT: Specific for the GST mu form (in rat and human) and does not cross-react with alpha or pi class GST's. Reactive with 0.1µg of purified rat and /or human GST Yb1.
- Application:** Western blot analysis.
- Quantity:** 100 µL of goat antiserum.
- Protein Concentration:** 68 mg/mL
- Working Dilution:**
A dilution of 1:1000 for Western blots will yield a positive result with > 0.1mg of purified rat Yb1 GST or recombinant human M1-1. Alkaline phosphatase-conjugated anti-goat IgG secondary antibody is used according to the manufacturer's instructions.
- Storage:** Stable for several weeks at 4°C, several months at -20°C, and indefinitely at -70°C. Avoid repeated freeze/thaw cycles.
- References:** Holmes, W. et.al., *Br. J. Cancer*, **62** : 209-212 (1990)