



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA
USA: 800-692-4633 • Fax: 248-852-4466
www.oxfordbiomed.com

Product specifications

**Preimmune Control Serum for
Polyclonal Anti-Human GST A1-1**

Product # GS 62c

Typical Lot

- Uses:** Use at same dilution as GS 62, (1:1000), for negative control with immunoassay or western blots
- Contents (per vial):** 100 μ L of preimmune rabbit serum.
- Stability:** Stable for several days at 4^o C, several months at -20^oC, and indefinitely at -70^oC. Avoid repeated freeze/thaw cycles.