



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA
USA: 800-692-4633 • Fax: 248-852-4466
www.oxfordbiomed.com

Product specifications

**Preimmune Control Serum for
Polyclonal Anti-Human GST P1-1**

Product # GS 72c Typical Lot

Applications:	Use at same dilution as GS 72, (1:1000), for negative control with immunoassay or western blots
Contents per Vial:	100 μ L of preimmune rabbit serum.
Stability:	Stable for several days at 4 ^o C, several months at -20 ^o C, and indefinitely at -70 ^o C. Avoid repeated freeze/thaw cycles.