

Product Specification Sheet

Description: Recombinant Human
15-Lipoxygenase Form-2
Product Number: LX26
Aliquot: 100 ng
Lot Number: lx26.050819CS
Storage: -70°C

DESCRIPTION:	Human recombinant 15-LOX-2 protein expressed from a baculoviral insect expression system.
MOLECULAR WEIGHT:	76 kDa
CONCENTRATION:	100 ng in 100 uL
STORAGE BUFFER:	10% SDS; 5% b-Mercatoethanol; Tris HCl; 0.01% Bromophenol Blue
STORAGE:	Stable for several months when stored at -20°C or for several years at -70°C.
PURITY:	Boil for 5 minutes and load 5-10 ng/ lane when using ECL or AP detection methods.
APPLICATIONS:	Western blotting standard.