



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA

USA 800-692-4633 • Fax 248-852-4466

www.oxfordbiomed.com

Product specifications

Native Ovine Prostaglandin H Synthase-1 Standard

Product # NP 01 Typical Lot

Contents: 10 µg of protein in 100 µL.

Source: Ovine seminal vesicles.

**Protein
Concentration:** 0.1 mg/mL.

Purity: Greater than 95% by electrophoresis.

Packaging: The enzyme is suspended in 100 mM Tris buffer (pH 8.0) containing 0.01% Tween.

Storage: Stable for several months at -20°C, and indefinitely at -70°C.

Recommended Use: Mix with sample buffer, boil, and run on SDS gel. For best results, use 1µg per lane when developing with alkaline phosphatase.