



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Human tPA Biotin Labeled IgG

Product Number: PA 33

Aliquot: 1.0 mg

Lot Number: pa33.1106A

Storage: -80°C

DESCRIPTION:	This rabbit polyclonal antibody is specific for human tissue plasminogen activator (tPA). The IgG preparation is provided with a biotin label.
MOLECULAR WEIGHT:	70 kDa
CONCENTRATION:	12.0 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.5 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
STORAGE:	-80°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Recombinant Human tPA
CLASS:	IgG
APPLICATIONS:	This antibody is suitable for use in Western blotting and ELISA application. Western: 1:300 dilution ELISA: 1:1000 dilution
SPECIES SPECIFICITY:	Human
REFERENCES:	Spence, MJ <i>et al.</i> <i>Cytokine</i> 18:1 pp 26-34 (April 7, 2002)