



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Human tPA Fluorescent Labeled IgG

Product Number: PA 34

Aliquot: 1.0 mg

Lot Number: pa34.804A

Storage: -70°C

DESCRIPTION:	This rabbit polyclonal antibody is specific for human tissue plasminogen activator (tPA). The IgG preparation has been fluorescent labeled to increase the utility of the antibody.
MOLECULAR WEIGHT:	70 kDa
CONCENTRATION:	1.7 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.5 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
STORAGE:	-80°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Recombinant Human tPA
CLASS:	IgG
APPLICATIONS:	This antibody is suitable for use in Western blotting and ELISA application. Western: 1:300 dilution ELISA: 1:1000 dilution
SPECIES SPECIFICITY:	Human