



OXFORD BIOMEDICAL RESEARCH
P.O. Box 522, Oxford MI 48371 • USA
USA: 800-692-4633 • Fax: 248-852-4466
www.oxfordbiomed.com

Product Specification

Anti-Human tPA (High Titer)

Product #: PA 35

Lot #: typical

Physical Appearance: yellow frozen solution

Solubility: > 10 mg/mL

Concentration: 0.6 mg/mL

Applications: Suitable for use in Western and EIA applications.

Considerations:

Storage Conditions: Store frozen at -70°C – stability > 1yr

AVOID REPEATED FREEZE THAW CYCLES

Buffer: 0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6

Other: < 98% pure by SDS PAGE