

Product Specification Sheet

Sheep Anti-Human tPA IgG
Fraction
Product Number: PA37
Aliquot: 1.0 mg
Lot Number: pa37.213
Storage: -70°C

DESCRIPTION:	Sheep Anti-tPA (Human), Sheep Polyclonal Antibody, IgG Fraction
ANTIGEN:	Recombinant human, full length, wild type protein
MOLECULAR WEIGHT:	70 kDa
CONCENTRATION:	2.94 mg/mL
SIZE:	1.0 mg
VOLUME:	0.34 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Sheep serum
BUFFER:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
PURIFICATION:	IgG fraction purified by immobilized Protein G

For Research Purposes Only.