



OXFORD BIOMEDICAL RESEARCH

## Product Specification Sheet

Human tPA Depleted Plasma (Na Citrate)

Product Number: PA43

Aliquot: 10.0 mL

Lot Number: Typical Lot

Storage: -70°C

<b>DESCRIPTION:</b>	Prepared from frozen human plasma using immobilized anti-human tPA IgG.
<b>MOLECULAR WEIGHT:</b>	N/A
<b>CONCENTRATION:</b>	N/A
<b>VOLUME:</b>	10.0 mL
<b>STORAGE:</b>	-70°C
<b>FORM:</b>	Frozen liquid
<b>SOURCE:</b>	Human plasma
<b>BUFFER:</b>	N/A

For research purposes only.