



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Human tPA, Immobilized
Product Number: PA48
Aliquot: 1.0 mg
Lot Number: Typical Lot
Storage: 4°C

DESCRIPTION:	May be used to immunopurify monoclonal and polyclonal antibodies directed against human tPA. May be used repeatedly
MOLECULAR WEIGHT:	70,000
CONCENTRATION:	2.0 mg/mL coupled
SIZE:	1.0 mg
VOLUME:	2.0 mL of a 1:1 slurry
STORAGE:	4°C
FORM:	Resin
SOURCE:	CHO cells
BUFFER:	0.1M HEPES; 0.1M NaCl; 0.02% Sodium Azide; pH 7.4

For research purposes only.