



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Single Chain Human tPA
Product Number: PA 51
Aliquot: 100 µg
Lot Number: pa51.506A
Storage: -80°C

DESCRIPTION: This single chain human tissue plasminogen activator-inhibitor (tPA) is recombinant human produced Chinese hamster ovary cells.

MOLECULAR WEIGHT: 75 kDa

CONCENTRATION: 3.9 mg/mL determined spectrophotometrically.

STORAGE BUFFER: 0.4 M HEPES; 0.1 M NaCl; pH 7.4

STORAGE: -80°C. **AVOID REPEATED FREEZE THAWS.**