



Product Specification Sheet

Rabbit Anti-Human tPA IgG
Fraction
Product Number: PA59
Aliquot: 1.0 mg
Lot Number pa59. Typical Lot
Storage: -70°C

DESCRIPTION:	Rabbit anti-human tPA, IgG Fraction
CONCENTRATION:	10 mg/ml
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
STORAGE:	-70°C
IMMUNOGEN:	Recombinant full length wt-Human tPA
MOLECULAR WEIGHT:	70 kDa
PURIFICATION:	Protein A purification
SPECIES SPECIFICITY:	Human