



OXFORD BIOMEDICAL RESEARCH

## Product Specification Sheet

Anti-Mouse tPA MoAB  
Product Number: PA 62  
Aliquot: 1.0 mg  
Lot Number: pa62.1003A  
Storage: -80°C

<b>DESCRIPTION:</b>	This is a capture monoclonal antibody produced in a tPA knockout mouse.
<b>MOLECULAR WEIGHT:</b>	70 kDa
<b>CLONE:</b>	H27B6
<b>CONCENTRATION:</b>	1.1 mg/mL determined spectrophotometrically.
<b>STORAGE BUFFER:</b>	0.05 M Sodium Phosphate; 0.15 M NaCl; pH 7.4
<b>STORAGE:</b>	-80°C
<b>IMMUNOGEN:</b>	Recombinant Mouse tPA
<b>PURIFICATION:</b>	Purified by Protein A Sepharose.
<b>CLASS:</b>	IgG
<b>PURITY:</b>	≥ 98% assessed by visual inspection of a Coomassie <sup>®</sup> blue-stained SDS-PAGE gel.
<b>APPLICATIONS:</b>	Suitable for use in ELISA, Western blot and immunohistochemistry. Recommended Western concentration at 1-5 ug/mL.
<b>SPECIES SPECIFICITY:</b>	Mouse and rat.