




OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Mouse tPA MoAB
Product Number: PA 63
Aliquot: 1.0 mg
Lot Number pa63.704A
Storage: -80°C

DESCRIPTION:	This is a monoclonal antibody produced in a tissue plasminogen activator (tPA) knockout mouse. The antibody reacts with both free and complexed forms of tPA.
MOLECULAR WEIGHT:	70 kDa
CLONE:	H6C5
CONCENTRATION:	2.3 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.15 M NaCl; pH 7.4
STORAGE:	-80°C
IMMUNOGEN:	Recombinant Mouse tPA
CLASS:	IgG
PURIFICATION:	Purified by Protein A Sepharose.
PURITY:	≥ 98% assessed by visual inspection of a Coomassie [®] blue-stained SDS-PAGE gel.
APPLICATIONS:	Suitable for use in ELISA and Western blot.
SPECIES SPECIFICITY:	Mouse



H6C5

1. Recombinant mouse tPA.

2. Recombinant mouse tPA in complex with PAI-1.