



Oxford Biomedical Research
Superior Science, Reliable Results.

Product Specification Sheet

Anti-Human uPA IgG Fraction

Product Number: PA 66

Aliquot: 1.0 mg

Lot Number pa66.706A

Storage: -80°C

DESCRIPTION:	This polyclonal antibody is produced in rabbit and is specific for human uPA. The antibody reacts with both free and complexed forms of uPA.
MOLECULAR WEIGHT:	54 kDa High Molecular Weight form 31 kDa Low Molecular Weight form
CONCENTRATION:	14.4 mg/mL determined Spectrophotometrically.
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.15 M NaCl; pH 7.4
STORAGE:	-80°C
IMMUNOGEN:	Recombinant Human uPA
CLASS:	IgG
PURIFICATION:	Purified by protein A sepharose.
APPLICATIONS:	Suitable for use in ELISA and Western blot.
SPECIES SPECIFICITY:	Human