

Product Specification Sheet

Description: Rabbit anti-Human uPA, Biotin Labeled
Product Number: PA67
Aliquot: 1.0 mg
Lot Number: pa67.409
Storage: -70°C

DESCRIPTION:	Rabbit anti-Human uPA, IgG fraction that has been biotin labeled.
MOLECULAR WEIGHT:	54 kDa
CONCENTRATION:	8.2 mg/mL
STORAGE BUFFER:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
STORAGE:	Store at -70°C
IMMUNOGEN:	Recombinant human wt full length glycosylated protein
PURIFICATION:	IgG purification using immobilized Protein A
APPLICATIONS:	ELISA, IHC, WB
SPECIES SPECIFICITY:	Human