

Product Specification Sheet

Description: Fluorescein Labeled Rabbit Anti-Human uPA (IgG Fraction)
Product Number: PA68
Aliquot: 1.0 mg
Lot Number: pa68.205
Storage: -70 °C

DESCRIPTION:	Fluorescein Labeled Rabbit Anti-Human uPA (IgG Fraction)
MOLECULAR WEIGHT:	160 kDa
CONCENTRATION:	7.4 mg/mL
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
STORAGE:	-70°C
IMMUNOGEN:	Recombinant human wild type, full length, glycosylated protein.
PURIFICATION:	Purified by immobilized Protein A.
CLASS:	IgG
APPLICATIONS:	ELISA, WB, IHC
SPECIES SPECIFICITY:	Human