



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Recombinant Active Rat tPA
Product Number: PA78
Aliquot: 0.1 mg
Lot Number: Typical Lot
Storage: -70°C

DESCRIPTION:	Recombinant Rat tPA produced in insect cells, >85% single chain, >95% active
MOLECULAR WEIGHT:	70,000 kD
CONCENTRATION:	2.0 mg/mL
SIZE:	0.1 mg
VOLUME:	0.05 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Insect cell culture
STORAGE BUFFER:	0.4M HEPES; 0.1M NaCl; pH 7.4

For research purposes only.