

**Product Specification Sheet**

Description: Human LMW uPA,  
Product Number: PA85  
Aliquot: 0.1 mg  
Lot Number: pa85.typical lot  
Storage: -70°C

<b>DESCRIPTION:</b>	This Human LMW Urokinase is cleaved from the amino terminal at Lys 135 and retains the proteinase domain.
<b>MOLECULAR WEIGHT:</b>	33 kDa
<b>CONCENTRATION:</b>	2.5 mg/mL
<b>STORAGE BUFFER:</b>	0.05 M Sodium acetate; 0.1 M NaCl; 1 mM EDTA; pH5.0
<b>STORAGE:</b>	Store at -70°C
<b>SOURCE:</b>	Insect cell culture, stable clone.
<b>PURITY:</b>	>95% by SDS-PAGE Analysis