

**Product Specification Sheet**

Description: Human LMW  
Urokinase, FITC conjugate  
Product Number: PA87  
Aliquot: 0.1 mg  
Lot Number: pa87.typical lot  
Storage: -70°C

**DESCRIPTION:** This Human LMW Urokinase is cleaved from the amino terminal at Lys 135 and retains the proteinase domain. It has been conjugated to FITC.

**MOLECULAR WEIGHT:** 33 kDa

**CONCENTRATION:** 2.5 mg/mL

**STORAGE BUFFER:** 0.05 M Sodium Acetate; 0.1 M NaCl; 1 mM EDTA; pH 5.0

**STORAGE:** Store at -70°C

**SOURCE:** Recombinant human urokinase produced in a stable insect cell clone.

**PURITY:** >95% by SDS-PAGE Analysis