



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA  
USA: 800-692-4633 • Fax: 248-852-4466  
[www.oxfordbiomed.com](http://www.oxfordbiomed.com)

### *Product specifications*

## Preimmune Control Serum for Polyclonal Anti-PGHS-1

Product # PG 19c      Typical Lot

- Uses:** Use at the same dilution as PG19 for negative control with immunoassay or Western blots.
- Contents(per vial):** 1 mL of preimmune goat anti-serum.
- Stability:** Stable for several days at 4 °C, several months at -20°C, and indefinitely at -70°C. Avoid repeated freeze/thaw cycles.
- Protein Concentration:** Approximately 5.4 mg/mL.