



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA

USA: 800-692-4633 • Fax: 248-852-4466

www.oxfordbiomed.com

Product specifications

Preimmune Control Serum for Polyclonal Anti-PGHS-1

Product # PG 20c Typical Lot

- Uses:** Use at same dilution as PG 20 for negative control with immunoassay or western blots.
- Contents (per vial):** 1 mL of preimmune rabbit anti-serum.
- Stability:** Stable for several days at 4° C, several months at -20°C, and indefinitely at -70°C. Avoid repeated freeze/thaw cycles.
- Protein Concentration:** Approximately 1 mg/mL