



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA

USA 800-692-4633 • Fax 248-852-4466

www.oxfordbiomed.com

Product specifications

Preimmune Control Serum for *Polyclonal Anti-PGHS-2*

Product # PG 26c Typical Lot

Protein Concentration: Approximately 1.1 mg/mL

Contents/vial: 1 mL of preimmune rabbit anti-serum

Applications: Use at same dilution as PG 26 for negative control with immunoassay or western blots.

Stability: Stable for several days at 4 °C. several months at -20°C. and indefinitely at -70°C. Avoid repeated freeze/thaw cycles.