

Product Specification Sheet

Description: Anti-PGF2a
Product Number: PG32
Aliquot: 100uL
Lot Number: 061600DM
Storage: -70C

DESCRIPTION: Anti-PGF2a, Sheep polyclonal antibody

MOLECULAR WEIGHT: 100 uL whole serum

STORAGE BUFFER: Neat serum with 0.09% sodium azide

STORAGE: Stable for at least one year when stored at -70°C

**RECOMMENDED
DILUTION:** 1:10,000

CROSS REACTIVITY:

PGF _{2a}	100%
6-keto-PGF _{1a}	0.9%
Thromboxane B ₂	16%
13,14-duhydro-15-keto PGF _{2a}	0.06%
PGD ₂	0.2%
PGE ₁	0.1%
PGE ₂	0.3%
6-keto PGE ₁	0.03%
PGF _{1a}	21%
Arachidonic acid	0.001%

SPECIES SPECIFICITY: Not species specific