



OXFORD BIOMEDICAL RESEARCH
P.O. Box 522, Oxford MI 48371 • USA
USA: 800-692-4633 • Fax: 248-852-4466
www.oxfordbiomed.com

Product Specification

Human Plasminogen Activator Inhibitor I (active wild type)

Product #: PI 01

Lot #: HWT-501

Molecular Weight: 43 KDa

Physical Appearance: colorless frozen solution

Solubility: > 2 mg/mL and < 5 mg/mL

Spectrophotometric Data:

Ultraviolet: Absorbance (280nm) = 2.65 0.1% = 1.0

Concentration = 2.65 mg/mL

Active Concentration = 2.65 mg/mL

Kinetic Data:

Second order rate constants for inhibition of:

UPA = $6.3 \times 10^6 \text{ M}^{-1}\text{s}^{-1}$

tPA = $5.7 \times 10^6 \text{ M}^{-1}\text{s}^{-1}$

Considerations:

Storage Conditions: Store frozen at -70°C – stability > 1yr

Buffer Composition: 0.05 M Sodium Phosphate; 0.1 M NaCl;

1 mM EDTA; pH 6.6

Other: > 98% pure by SDS PAGE 95 (+/-) 5% active by uPA titration

