

Human PAI-1 Latent Wild Type

Product Number: PI02

Aliquot: 0.5 mg

Lot Number: pi02.805

Storage: -70°C

DESCRIPTION: This preparation is <0.4% active as determined by titration with HMW tc-urokinase and SDS PAGE. Latent PAI-1 is a useful control for experiments when used in conjunction with the active form above.

MOLECULAR WEIGHT: 43 kD

CONCENTRATION: 2 mg/mL


STORAGE BUFFER: 0.1 M HEPES; 0.1 M NaCl; pH 7.4

STORAGE: -70°C. Avoid repeated Freeze-Thaw cycles.

SOURCE: E. coli

FORM: Frozen liquid

PURITY: ≥98% by SDS-PAGE gel.



10% SDS-PAGE - GelCode Blue Stain

1 PAI-L-697 (5 ug) Reduced

2 PAI-L-697 (5 ug)+uPA-HTC (8 ug) Reduced

3 Prestained Standard

For in vitro research use only.