



Product Specification Sheet

Human PAI-1 (Stable Mutant)
Product Number: PI 05
Aliquot: 0.5 mg
Lot Number pi05.803
Storage: -80°C

DESCRIPTION: Human Plasminogen Activator Inhibitor 1 (Stable Mutant)

CONCENTRATION: 2.5 mg/mL

MOLECULAR WEIGHT: 43,000

STORAGE BUFFER: 0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6

STORAGE: Store at -80°C. Avoid multiple freeze/thaw cycles.

VOLUME: 0.2 mL