

Product Specification Sheet

Description: Rat PAI-1
Product Number: PI20
Aliquot: 0.5 mg
Lot Number: pi20.1099
Storage: -70°C

DESCRIPTION: Wild-type latent Rat PAI-1 produced in *E. coli*

MOLECULAR WEIGHT: 43 kDa

CONCENTRATION: 2.1 mg/mL

STORAGE BUFFER: 0.1 M HEPES; 0.1 M NaCl; pH 7.4

STORAGE: -70°C

PURITY: >95% by SDS-PAGE analysis