



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-PAI-1 (Human) Rabbit Polyclonal
Antisera
Product Number: PI30
Aliquot: 1.0 mL
Lot Number: 1002
Storage: -70°C

DESCRIPTION:	Anti-PAI-1 (Human) Rabbit Polyclonal Antisera
STORAGE:	-70°C
IMMUNOGEN:	Human PAI-1
PURIFICATION:	None
VOLUME:	1.0 mL
APPLICATIONS:	ELISA, IHC, WB
SPECIES SPECIFICITY:	Human